



**Panc02-Luc | 305706**

**Protein expression**

**XXXXXXXXXX**

**Culture Medium** RPMI 1640, w: 2.0 mM , w: 2.0 /  NaHCO3 (  Cytion 820700a)

**Supplements**  FBS 10%

**Dissociation Reagent**

**Doubling time** 24-48

**Subculturing**  PBS  T25

**Seeding density** 1  3 <sup>4</sup> /

**Fluid renewal** 2  3

**Freeze medium**  + 10% DMSO

- Thawing and Culturing Cells**
- - cryovial  -150
  - 37
  - 70%
  - 15
  - 200 x g  5
  -

