

NCI-H820 | 305841

Thawing and Culturing Cells

1. **Thawing**
2. **Resuspension** cryovial **-150**
3. **Incubation** **37°C**
4. **Seeding** **70%**
5. **Medium** **15**
6. **Centrifugation** **300 x g** **3**
7. **Resuspension** **10**
8. **Seeding**

Incubation Atmosphere 37°C, 5% CO₂,

Flask Coating

Shipping Conditions

Storage Conditions **-150** **-196 °C**

/ / HLA

Sterility **PCR**