

Product sheet

KHYG-1 | 305890

Mutational profile p.Gly12Ala, ; p.Arg248Trp,

Culture Medium RPMI 1640, w: 2.0 mM , w: 2.0 / NaHCO₃ (Cytion 820700a)

Supplements FBS 10% IL-2 10

Dissociation Reagent

Doubling time 24-48 ; ~30-40 ; ~54 , ~30 , ~25

Split ratio 1/4 3-4

Fluid renewal 3-4

Freeze medium + 10% DMSO

- Thawing and Culturing Cells**
- 1.
 2. cryovial -150
 3. 37
 4. 70%
 5. 15
 6. 200 x g 5
 - 7.

Incubation Atmosphere 37°C, 5% CO₂,

