

CHO-CD36 | 305979

CHO-CD36

Description

CHO-CD36 (CHO) CD36 (oxLDL), -1, CD36

CHO-CD36

Organism CHO

Tissue CHO

CHO

Age CHO

Gender CHO

Morphology CHO

Cell type CHO

CHO-CD36

Citation CHO-CD36 (Cytion 305979)

Biosafety level 1

NCBI_TaxID 10029

CellosaurusAccession CVCL_8848

CHO-CD36

Receptors expressed CD36

CHO-CD36 | 305979

CHO-CD36

Culture Medium DMEM:Ham's F12 (1:1), w: 3.1 g/l / w: 2.5 g/l L-Asparagine, w: 15 g/l
 CHO A (InSCREENeX; InSCREENeX I)

Supplements FBS 5% Geneticin (G418-Sulfat) 0.5 mg/ml

Dissociation Reagent EDTA

Subculturing PBS

Fluid renewal 2-3

Post-Thaw Recovery 1:2 1:3 T25

Freeze medium (FBS) + 10% DMSO

- Thawing and Culturing Cells**
1. Thaw cryovial in 37°C water bath
 2. Centrifuge cryovial -150 rpm
 3. Wash cells 37°C
 4. Resuspend cells 70%
 5. Seed cells 15
 6. Spin cells 300 x g 3 min
 7. Wash cells 10
 8. Resuspend cells

