

CHO-NECTIN4 | 305984

CHO-NECTIN4

Description	CHO-NECTIN4 (CHO) Nectin-4, ADCs,
Organism	
Tissue	
Disease	NECTIN4 (PVRL4)
Applications	ADC; NECTIN4;

CHO-NECTIN4

Age	
Gender	
Morphology	
Cell type	
Growth properties	

CHO-NECTIN4

Citation	CHO-NECTIN4 (Cytion 305984)
Biosafety level	1
NCBI_TaxID	10029
CellosaurusAccession	CVCL_A8W9

CHO-NECTIN4 | 305984

GMO Status

GMO-S1: CHO NECTIN4

Surface antigens

NECTIN4 (NECTIN4/112)

Culture Medium

DMEM:Ham's F12 (1:1), w: 3.1 / w: 2.5 L- , w: 15
CHO A (InSCREENeX; InSCREENeX I

Supplements

FBS 5% Geneticin (G418-Sulfat) 0.5

Dissociation Reagent

EDTA

Doubling time

14-16

Subculturing

1:2 1:3 T25 PB

Split ratio

1:5

Seeding density

2×10^4 /

Fluid renewal

2-3

Post-Thaw Recovery

1:2 1:3 T25

Freeze medium

(FBS) + 10% DMSO

CHO-NECTIN4 | 305984

**Thawing and
Culturing Cells**

1. Thaw the cells quickly in a water bath at 37°C.
2. Add the cells to a cryovial and store at -150 °C.
3. Thaw the cells quickly in a water bath at 37°C.
4. Add the cells to a cryovial and store at -150 °C.
5. Add the cells to a cryovial and store at -150 °C.
6. Add the cells to a cryovial and store at -150 °C.
7. Add the cells to a cryovial and store at -150 °C.
8. Add the cells to a cryovial and store at -150 °C.

**Incubation
Atmosphere**

37°C, 5% CO₂, humidified

**Shipping
Conditions**

Cells are shipped in a cryovial with dry ice.

**Storage
Conditions**

Cells are stored at -150 °C to -196 °C.

Cells are not for use in clinical trials.

Sterility

Cells are tested for sterility by PCR.

Cells are tested for sterility by PCR.