

| 300695

Description	CTIGM.Melano Growth Medium for Melanocytes		
Organism			
Tissue			
Disease			
Metastatic site			
Applications			
Age	1-17		
Gender			
Ethnicity			
Morphology			
Cell type			
Growth properties			
Citation			
Biosafety level	1		
NCBI_TaxID	9606		

| 300695

CellosaurusAccession XXXXXXXXXX

GMO Status XXXXXXXXXXXXXXXXXXXXXXXXXXXX; XXXXXXXXXXXXXXXXXXXXXXX

Culture Medium CTIGM.Melano : XXXXXXXXXXXXXXXXXXXX

Supplements XXXXXXXXXXXX CTIGM XXXXX XXXXXXXXXXXXXXXXXXXXXXXXXXXX Melano Growth Medium

Dissociation Reagent XXXXXXXXXX

Doubling time XXXXXX 36 XXX 48 XXXXXXXXXX

Subculturing XXXXXXXXXXXXXXXXXXXX XXXXXXXXXX PBS XXXXXXXXXXXXXXXXXXXXXXXXXXXX XXXXXXXXXX Accutase XXXXXXXXXXXXXXXXXXXX 37°C XXXXXXXXXX

Split ratio 1 XXX 3

Seeding density 1 XXX 2×10^3 XXXXXX/XX²

Fluid renewal XXX 2 XXX 3 XXX

Freeze medium XXXXXXXXXXXXXXXXXXXXXXXXXXXX XXXXXXXXXXXXXXXXXXXXXXXXXXXX (XXXXXX FBS) + 10% DMSO

| 300695

Thawing and
Culturing Cells

- 1. Thaw cells rapidly in a water bath at 37°C.
- 2. Mix the cells gently and transfer to a cryovial.
- 3. Store the cells at -150 °C.
- 4. Thaw cells rapidly in a water bath at 37°C.
- 5. Mix the cells gently and transfer to a cryovial.
- 6. Store the cells at -150 °C.
- 7. Thaw cells rapidly in a water bath at 37°C.
- 8. Mix the cells gently and transfer to a cryovial.

Incubation
Atmosphere

37°C, 5% CO₂, humidified

Shipping
Conditions

Cells are shipped on dry ice at -196 °C.

Storage
Conditions

Cells are stored at -150 °C.