

DAKIKI | 305928

Description	DAKIKI (EBV) IgA1 IgA1 DAKIKI Golgi DAKIKI CMP-NeuAc:GalNAc-IgA1 α 2,6-sialyltransferase GalNAc IgA1 DAKIKI
Organism	
Tissue	
Synonyms	, 1
Age	
Gender	
Ethnicity	
Morphology	
Cell type	
Growth properties	
Citation	(Cytion 305928)
Biosafety level	2
NCBI_TaxID	9606
CellosaurusAccession	CVCL_3675

DAKIKI | 305928

Culture Medium	RPMI 1640, w: 2.0 mM , w: 2.0 / NaHCO3 (Cytion 820700a)
Supplements	 FBS 10%
Dissociation Reagent	
Subculturing	 ⁵
Seeding density	2 ⁵ ²
Fluid renewal	2 3
Freeze medium	 + 10% DMSO
Thawing and Culturing Cells	1. 2. cryovial -150 3. 37 4. 70% 5. 15 6. 200 x g 5 7.
Incubation Atmosphere	37°C, 5% CO ₂ ,

DAKIKI | 305928

Shipping Conditions

Storage Conditions

-150 -196 °C