



XXXXXXXXXXXXXXXXXXXX





XXXXXX

XXXXXXXXXXXXXXXXXXXX

C57BL/6


XXXXXX

GL-261-Luc (XXXXXXXXXXXXXXXXXXXX Cytion 305662)

1

10090

CVCL_C9CB

GMO-S1:  (Luc2, ) 

333

Luc2 (XXXXXXXX, XXXXXX-XXXXXXXXXXXX)

Product sheet

GL261-Luc | 305662

General information

Culture Medium DMEM, w: 4.5 g/l Glucose, w: 4 g/l Sodium butyrate, w: 3.7 g/l NaHCO₃, w: 1.0 g/l Penicillin, w: 1.0 g/l Streptomycin, w: 1.0 g/l Fungizone

Supplements FBS 10%

Dissociation Reagent Trypsin-EDTA

Subculturing PBS, 10 min, 37°C, T25 flasks

Seeding density 1 x 10⁴ cells / cm²

Fluid renewal 2 x 3 days

Freeze medium FBS + 10% DMSO

Thawing and Culturing Cells

1. Thaw cells quickly in a 37°C water bath.
2. Add cells to a cryovial and store at -150°C.
3. Thaw cells quickly in a 37°C water bath.
4. Add cells to a flask and incubate at 37°C, 5% CO₂.
5. Add cells to a flask and incubate at 37°C, 5% CO₂.
6. Add cells to a flask and incubate at 37°C, 5% CO₂.
7. Add cells to a flask and incubate at 37°C, 5% CO₂.

Incubation Atmosphere 37°C, 5% CO₂

Shipping Conditions 4°C

GL261-Luc | 305662

—150 —196 °C

