

NOTE: STI STI

A375-GFP A375



Western blot analysis showing eGFP and A375-GFP expression in vitro. The blot displays two rows of bands. The top row is labeled 'eGFP' and the bottom row is labeled 'A375-GFP'. There are four lanes in total, with the first two lanes corresponding to the eGFP row and the next two to the A375-GFP row. The bands in the A375-GFP row are significantly more intense than those in the eGFP row, indicating higher expression levels of A375-GFP.

A375-GFP

XXXXXX

XX, XXXXXXXX

XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX

54 田田

XXXXXXXXXX

XXXXXXXXXX

A375-GFP (Cytion 305665)

1

9606

CVCL_QZ67

GM0-S1: ZsGreen1

Cytion GmbH | XXXXXXXX-XXXX3 | 69124 XXXXXXXX | XXXXXXXX
 XXXX +49(0)6221 405780 | www.cytion.com | info@cytion.com

Product sheet

A375-GFP | 305665


Incubation
Atmosphere

37°C, 5% CO₂, 

Shipping
Conditions



Storage
Conditions

 -150  -196 °C 