

HEK293-Rpn11-HTBH | 305719

HEK293-Rpn11-HTBH

Description	HEK293 - Rpn11-HTBH HEK293 JAMM (interactome studies) 26S (UPS)
Organism	
Tissue	
Disease	/ (HEK293; Rpn11-HTBH)
Applications	26S; Rpn11/PSMD14; -; (hexahistid

HEK293

Morphology	
Cell type	
Growth properties	

HEK293-Rpn11-HTBH

Citation	HEK293 -Rpn11-HTBH (Cytion 305719)
Biosafety level	1
NCBI_TaxID	9606
GMO Status	GMO-S1: HEK293 -Rpn11-HTBH (hexahistid

HEK293-Rpn11-HTBH

HEK293

Product sheet

HEK293-Rpn11-HTBH | 305719

Culture Medium	DMEM, w: 4.5 g/l Glucose, w: 4 g/l L-Glutamine, w: 3.7 g/l NaHCO ₃ , w: 1.0 g/l Penicillin, w: 1.0 g/l Streptomycin
Supplements	10% FBS
Dissociation Reagent	PBS + Accutase 5 min
Doubling time	24-36 h
Subculturing	PBS + Accutase
Split ratio	1:10
Seeding density	2 × 10 ⁴ cells/cm ²
Fluid renewal	2-3 times
Freeze medium	DMEM + 10% DMSO
Thawing and Culturing Cells	1. Thaw cells in a 37°C water bath. 2. Add cells to a cryovial and store at -150 °C. 3. Add cells to a 37 °C water bath. 4. Add cells to a 70% ethanol solution. 5. Add cells to a 15% ethanol solution. 6. Add cells to a 200 x g centrifuge tube. 7. Add cells to a 5 ml tube.
Incubation Atmosphere	37°C, 5% CO ₂ , humidified

Product sheet

HEK293-Rpn11-HTBH | 305719

Shipping
Conditions

Storage
Conditions

-150 -196 °C