

BC-3C | 305246

BC-3C

Description	BC3c 12 H-ras 12 13
Organism	
Tissue	
Disease	
Synonyms	BC3c

BC-3C

Age	82
Gender	
Ethnicity	

Growth properties

BC-3C

Citation	BC-3C (Cytion 305246)
Biosafety level	1
NCBI_TaxID	9606
CellosaurusAccession	CVCL_1958

BC-3C

BC-3C

Product sheet

BC-3C | 305246

Culture Medium McCoy's 5a, w: 3.0 g/l / 3.0 g/l, w: 2.0 g/l / 2.0 g/l, w: 2.2 g/l / 2.2 g/l NaHCO₃ (g/l)

Supplements FBS 10%

Dissociation Reagent Accutase 20 mg/ml 37°C;

Doubling time ~25-30 days

Freeze medium + 10% DMSO

Thawing and Culturing Cells

1. Thaw cells in a 37°C water bath.
2. Centrifuge cells at 300 x g for 5 min.
3. Wash cells with PBS.
4. Resuspend cells in 1 ml of culture medium.
5. Seed cells into a 25 cm² flask.
6. Incubate cells at 37°C, 5% CO₂.
7. Monitor cell growth.

Incubation Atmosphere 37°C, 5% CO₂

Shipping Conditions

Storage Conditions -150 °C to -196 °C