

JHH-7 | 305386

General information

Description	JHH-7 (HCC) University) HCC JHH-7 (HCC) HCC, TP53
Organism	
Tissue	
Disease	
Metastatic site	
Applications	
Synonyms	Jhh-7, JHH7, FLC-7, FLC7, -7

Characteristics

Age	53
Gender	
Ethnicity	
Morphology	
Cell type	
Growth properties	

Additional information

Citation	JHH-7 (Cytion 305386)
Biosafety level	1

XXXXXXXX JHH-7 | 305386

NCBI_TaxID 9606

CellosaurusAccession CVCL_2805

XXXXXXXXXXXXXXXXXXXX

Mutational profile XXXXXXXXXXXXXXXX p.Arg267Pro, XXXXXXXX

XXXXXXXXXXXX

Culture Medium Williams' E medium XXXX MEM + 10% FBS (XX)

Supplements XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX FBS 10%

Dissociation Reagent 0.25% Trypsin-EDTA XXXX Accutase

Doubling time 26.4 XXXXXXXX; 17.1 XXXXXXXX; 30 XXXXXXXX, ~29 XXXXXXXX, ~26 XXXXXXXX

Subculturing XXXXXXXXXXXX XXXXXXXXXXXX PBS XX XXXXXXXXXXXX Trypsin-EDTA 0.25% XXXX Accutase XXXXXXXXXXXXXXXXXXXXXXXX

Split ratio 1 XXX 5

Seeding density 1 XXX 3 × 10⁴ XXXXX/XXXXXXXXXXXXXXXXXXXX

Fluid renewal XXX 2 XXX 3 XXX

Freeze medium XX XX + 10% DMSO XXXXXXXXXXXXXXXX

JHH-7 | 305386

**Thawing and
Culturing Cells**

1.

Thaw the cells in a water bath at 37°C. Transfer the cells to a cryovial and store at -150 °C.
2.

Remove the cells from the cryovial and transfer to a tube. Centrifuge at 200 x g for 5 min.
3.

Resuspend the cells in 15 ml of medium. Seed the cells into a 70% confluent well.
4.

Incubate the cells at 37°C, 5% CO₂.
5.

After 15 days, the cells should be ready for use.
6.

After 200 x g for 5 min, the cells should be ready for use.
7.

After 15 days, the cells should be ready for use.

**Incubation
Atmosphere**

37°C, 5% CO₂, **JHH-7**

**Shipping
Conditions**

JHH-7

**Storage
Conditions**

JHH-7 -150 °C -196 °C **JHH-7**