

Product sheet

T Hs 611 | 305915

|                   |                                                                                         |
|-------------------|-----------------------------------------------------------------------------------------|
| Description       | Hs 611.T <span></span> (cHL) <span></span><br>Hs 611.T <span></span> HRS, <span></span> |
| Organism          | <span></span>                                                                           |
| Tissue            | <span></span>                                                                           |
| Disease           | <span></span>                                                                           |
| Metastatic site   | <span></span>                                                                           |
| Applications      | <span></span> ; <span></span> NF-κB <span></span> JAK-STAT; <span></span>               |
| Synonyms          | Hs-611-T, Hs611T, HS611T                                                                |
| Age               | 48 <span></span>                                                                        |
| Gender            | <span></span>                                                                           |
| Ethnicity         | <span></span>                                                                           |
| Morphology        | <span></span>                                                                           |
| Cell type         | <span></span>                                                                           |
| Growth properties | <span></span>                                                                           |
| Citation          | Hs 611.T ( <span></span> Cytion 305915)                                                 |
| Biosafety level   | 1                                                                                       |
| NCBI_TaxID        | 9606                                                                                    |

T Hs 611 | 305915

CellosaurusAccession CVCL\_0823

**Culture Medium** DMEM, w: 4.5 g/l Glucose, w: 4 g/l Sodium Pyruvate, w: 3.7 g/l NaHCO<sub>3</sub>, w: 1.0 g/l Penicillin, w: 1.0 g/l Streptomycin

**Supplements** FBS 10%

**Dissociation Reagent** Trypsin

**Doubling time** 24 h 48 h

**Subculturing** 3–5 × 10<sup>5</sup> cells / 100 cm<sup>2</sup> flask

**Split ratio** 1 : 3

**Seeding density** 2 × 10<sup>4</sup> cells / 100 cm<sup>2</sup> flask

**Fluid renewal** 2–3 times

**Freeze medium** DMEM + 10% DMSO

T Hs 611 | 305915

Thawing and  
Culturing Cells

- 1. Thaw the cells in a water bath at 37°C.
- 2. Add the cells to a cryovial and store at -150 °C.
- 3. Thaw the cells in a water bath at 37°C.
- 4. Add the cells to a cryovial and store at 70%.
- 5. Add the cells to a cryovial and store at 15 °C.
- 6. Add the cells to a cryovial and store at 200 x g 5 min.
- 7. Add the cells to a cryovial and store at 15 °C.

Incubation  
Atmosphere

37°C, 5% CO<sub>2</sub>, humidified

Shipping  
Conditions

Store at -150 °C.

Storage  
Conditions

Store at -150 °C to -196 °C.