

CHO-K1 () | 603486

Description	CHO-K1 CHO 1950 CHO-K1 CHO-K1 CHO-K1 CHO-K1 CHO
-------------	---

Organism	
----------	--

Tissue	
--------	--

Applications	
--------------	--

Synonyms	CHO K1, CHOK1, CHO cell clone K1, GM15452
----------	---

Age	
-----	--

Gender	
--------	--

Morphology	
------------	--

Growth properties	
-------------------	--

Citation	CHO-K1 (Cytion 603486)
----------	------------------------

Biosafety level	1
-----------------	---

NCBI_TaxID	10029
------------	-------

CellosaurusAccession	CVCL_0214
----------------------	-----------

Product sheet


CHO-K1 () | 603486

Virus susceptibility	(), 2, ,
Reverse transcriptase	
Karyotype	50 2n = 22.

Culture Medium	CD CHO ThermoFisher,
Supplements	L- 8 mM (40 mL 1 200 mM)
Dissociation Reagent	
Doubling time	22
Subculturing	5/
Seeding density	4 × 10 ⁵ /
Fluid renewal	2 3
Freeze medium	(FBS) + 10% DMSO

XXXXXXXX CHO-K1 (XXXXXXXXXXXX) | 603486

Thawing and Culturing Cells

1.  
2.   cryovial  -150  
3.   37°C 
4.   70%
5.   15 
6.  300 x g  3  
7.   10  
8. 

Incubation Atmosphere

37°C, CO₂ 5%,    Erlenmeyer 

Shipping Conditions

Storage Conditions

 -150 -196 °C

[illegible]

Sterility