

A704 | 300217

Thawing and Culturing Cells

- 1.
2. cryovial -150
3. 37°C
4. 70%
5. 15
6. 300 x g 3
7. 10
- 8.

Incubation Atmosphere

37°C, 5% CO₂,

Flask Coating

Freezing Procedure

Shipping Conditions

Storage Conditions

-150 -196 °C

/ / HLA

Sterility

PCR

A704  | 300217

 HLA **A***: '34:02:01, '74:01:01
 B*: '35:01:01, '44:03:01
 C*: '04:01:01
 DRB1*: 15:03:01G
 DQA1*: '01:02:01
 DQB1*: '06:02:01
 DPB1*: '02:01:19, '04:02:01G
 E: '01:01:01, '01:03