

XXXXXXXXXXXXXXXXXXXX

XXXXXX

XX

XXXXXXXXXXXXXXXXXXXXXXXXXXXX

Product sheet

XXXXXXWiDr | 300377

Antigen expression	HLA A24, A32, B15, B18
Isoenzymes	PGM1, 1-2, PGM3, 1-2, G6PD, B, ES-D, 1, PEP-D, 1, 6PGD, A
Tumorigenic	XXX XXXXXXXXXXXXXXXXXX
Products	XXXXXXXXXXXXXXXXXXXXXXXXXXXX (CEA) 118 XXXXXXXX/10^6 XXXXX/10 XXX, XXXXXXXXXXXXXXXXXXXX (CSAp), XXXXXXXXXXXXXXXXXX
XXXXXXXXXXXXXX	
Culture Medium	DMEM:Ham's F12 (1:1), w: 3.1 XXXXX/XXXX XXXXXXX, w: 2.5 XXXXXXXXXXXXX L-XXXXXXX, w: 15 XXXXXXXXXXXXX HEPES, w: 0.5 XXXXXXXXXXXX
Supplements	XXXXXXXXXXXXXXXXXXXXXXXXXXXX FBS 10%
Dissociation Reagent	XXXXXXXXXX
Subculturing	XXX PBS XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX XXXXXXXXXXXX T25 XXXXXXX
Seeding density	1 x 10 ⁴ XXXXX ² /XX [^]
Fluid renewal	1 XXXX 2 XXXXXXXXXXXXXXXXXX
Post-Thaw Recovery	XXXXXXXXXX ⁴ XXXXXXXXXXXXXXXXXX XXX 24 XXXXXXX
Freeze medium	XX XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX (XXXXXXX FBS) + 10% DMSO

Thawing and Culturing Cells

1.              

Incubation Atmosphere

37°C, 5% CO₂, 

Flask Coating

☐ ☐ ☐ ☐ ☐

Freezing Procedure

Shipping Conditions

Storage Conditions

-150 -196 °C

Sterility

PCR

[illegible]

XXXXXXXXWiDr | 300377

XXXXXXXXXXXXHLA **A***: '01:01:01, '24:03:01
 B*: '35:01:01, '44:03:01
 C*: '04:01:01
 DRB1*: '04:02:01, '07:01:01
 DQA1*: '02:01:01, '03:01:01
 DQB1*: '02:02:01, '03:02:01
 DPB1*: '04:01:01
 E: '01:01:01, '01:03:01