

XXXXXXXXXX H-4-II-E | 305204

XXXXXXXXXXXXXXXXXXXX

Description	XXXXXXXXXX H-4-II-E XX (hepatoma) XXXXXXXX XX XXXXXXXXXXXXXXXXXXXX XXXXX H-4-II-E XX XX
Organism	XXX
Tissue	XXX
Disease	XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
Synonyms	H-4-II-E, H4IIE

XXXXXXXXXX

Breed/Subspecies	XXXXXXXXXX
Gender	XXX
Morphology	XXXXXXXXXXXXXXXXXXXX
Growth properties	XXXXXXXXXX

XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX

Citation	H-4-II-E (XXXXXXXXXXXXXXXXXXXXXXXXXXXX Cytion 305204)
Biosafety level	1
NCBI_TaxID	10116
CellosaurusAccession	CVCL_0284

XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX

XXXXXXXXXXXXXX

Product sheet

HEK293T H-4-II-E | 305204

Culture Medium	EMEM (MEM Eagle), w: 2 mM L-Glutamine, w: 2.2 g/L NaHCO ₃ , w: EBSS (Cytion 820100a)
----------------	---

Supplements FBS 10% NEAA 1%

Dissociation Reagent	
----------------------	--

Subculturing PBS T25

Seeding density	2 × 10 ⁴ cells/cm ²
-----------------	---

Fluid renewal 2 × 3

Freeze medium	(FBS) + 10% DMSO
---------------	------------------

Thawing and Culturing Cells

- 1.
2. cryovial -150
3. 37°C
4. 70%
5. 15
6. 300 x g 3
7. 10
- 8.

Incubation Atmosphere	37°C, 5% CO ₂ ,
-----------------------	----------------------------

Shipping Conditions

XXXXXX H-4-II-E | 305204

Storage Conditions

—150 —196 °C

Sterility

The diagram illustrates a PCR assay. It features a long horizontal bar with a repeating cross-hatch pattern. In the center of this bar, the letters "PCR" are displayed in a bold, black, sans-serif font. The bar is divided into three sections by the "PCR" label: a long left section, a short middle section containing the label, and a long right section. The entire bar is set against a light gray background.

[illegible]