

XXXXXXXXXXXXXXXXXXXX

XXXXXXXXXX

XX

XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX

Product sheet

HEK293T M2-10B4 | 400428

Products	HEK293T, HEK293T M2-10B4 4 (HEK293T M2-10B4 1 (-) HEK293T M2-10B4 8 (-))
HEK293T M2-10B4	
Culture Medium	RPMI 1640, w: 2.0 mM β -mercaptoethanol, w: 2.0 mM / 10 mM NaHCO ₃ (HEK293T M2-10B4 Cytion 820700a)
Supplements	HEK293T M2-10B4 FBS 10%
Dissociation Reagent	HEK293T M2-10B4
Subculturing	HEK293T M2-10B4 PBS HEK293T M2-10B4 T25 HEK293T M2-10B4
Seeding density	1×10^4 HEK293T M2-10B4 / well
Fluid renewal	2 HEK293T M2-10B4 HEK293T M2-10B4
Post-Thaw Recovery	HEK293T M2-10B4
Freeze medium	HEK293T M2-10B4 HEK293T M2-10B4 (HEK293T M2-10B4 FBS) + 10% DMSO

XXXXXX M2-10B4 | 400428

Thawing and Culturing Cells

1.              

Incubation Atmosphere

37°C, 5% CO₂, 

Flask Coating

XXXXXX

Freezing Procedure

Shipping Conditions

Storage Conditions

-150 -196 °C

Sterility

PCR