

MeWo | 300285

Thawing and Culturing Cells

1. **Thawing** **MeWo** cryovial **MeWo** -150 **MeWo** **MeWo**
2. **MeWo** cryovial **MeWo** -150 **MeWo** **MeWo**
3. **MeWo** **MeWo** 37°C **MeWo**
4. **MeWo** **MeWo** 70%
5. **MeWo** **MeWo** 15 **MeWo**
6. **MeWo** 300 x g **MeWo** 3 **MeWo** **MeWo**
7. **MeWo** **MeWo** 10 **MeWo** **MeWo**
8. **MeWo**

Incubation Atmosphere

37°C, 5% CO₂, **MeWo**

Flask Coating

MeWo **MeWo**

Freezing Procedure

MeWo

Shipping Conditions

MeWo

Storage Conditions

MeWo -150 **MeWo** -196 °C **MeWo**

MeWo / **MeWo** / HLA

Sterility

MeWo PCR **MeWo**

MeWo **MeWo** **MeWo**

XXXXXXXX MeWo | 300285

XXXXXXXXXXXX HLA **A***: '02:01:01, '26:01:01
 B*: 14:02:01, 38:01:01
 C*: '08:02:01, '12:03:01
 DRB1*: '01:02:01, '11:01:01G
 DQA1*: '01:01:02, '05:05:01
 DQB1*: '03:01:01G, '05:01:01G
 DPB1*: '04:01:01G, '04:02:01G
 E: '01:xx, '01:03:01