

SK-MEL-5 | 300157

Thawing and Culturing Cells

1. [REDACTED]
2. [REDACTED] cryovial [REDACTED]-150 [REDACTED]
3. [REDACTED] 37°C [REDACTED]
4. [REDACTED] 70%
5. [REDACTED] 15 [REDACTED]
6. [REDACTED] 300 x g [REDACTED] 3 [REDACTED]
7. [REDACTED] 10 [REDACTED]
8. [REDACTED]

Incubation Atmosphere 37°C, 5% CO₂, [REDACTED]

Flask Coating [REDACTED]

Freezing Procedure [REDACTED]

Shipping Conditions [REDACTED]

Storage Conditions [REDACTED]-150 [REDACTED]-196 °C [REDACTED]

[REDACTED] / [REDACTED] / HLA

Sterility [REDACTED] PCR [REDACTED]
[REDACTED]

XXXXXXXXXX SK-MEL-5 | 300157

XXXXXXXXXX STR **Amelogenin:** XXXXXXXX, XXXXXXXX
CSF1PO: 10,13
D13S317: 10,12
D16S539: 10,12
D5S818: 11,13
D7S820: 9,12
TH01: 6,9
TPOX: 11
vWA: 14,18
D3S1358: 16,17
D21S11: 29
D18S51: 15,16
Penta E: 5,12
Penta D: 9,11
D8S1179: 12,15
FGA: 20,2,22

XXXXXXXXXXXX HLA **A*:** '02:01:01, '11:01:01
B*: '07:02:01, '40:01:02
C*: '03:04:01, '07:02:01
DRB1*: '04:01:01, '13:01:01
DQA1*: '01:03:01, '03:01:01
DQB1*: '03:02:01, '06:03:01
DPB1*: '03:01:01, '16:01:01
E: '01:01, '01:03