

XXXXXXXXXXXXXXXXXXXX

XXXXXXXXXX

[illegible][illegible]


Product sheet

XXXXXX KG-1a | 300234

Antigen expression	HLA A30, A31, B35, Cw4
Isoenzymes	G6PD, B, PGM1, 1-2, PGM3, 0, ES-D, 1, Me-2, 1, AK-1, 0, GLO-1, 2
Viruses	EBNA (EBNA): XX
Reverse transcriptase	XXXXXX
XXXXXXXXXXXX	
Culture Medium	IMDM, w: 4.5 XXXX/XXXX XXXXX, w: 4 XXXXXXXXXXXX L-XXXXXXX, w: 25 XXXXXXXXXXXX HEPES, w: 1.0 XXXXXXXXXXXX XXXXXXXXXXXX
Supplements	XXXXXXXXXXXXXXXXXXXXXXX FBS 10%
Doubling time	45 XXXXXXX
Subculturing	XXXXXXXXXXXXXXXXXXXXXXX 5 XXXXXXX, 6 XXXXXXX, XXXXXXXXXXXX 300xg XXXXXXXX 3 XXXXXXX
Fluid renewal	XXX 3 XXX
Post-Thaw Recovery	XX 24 XXXXXXX
Freeze medium	XX XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX (XXXXXX FBS) + 10% DMSO

XXXXXX KG-1a | 300234

Thawing and Culturing Cells

1.  
2.   cryovial  -150  
3.   37°C 
4.   70%
5.   15 
6.  300 x g  3  
7.   10   
8. 

Incubation Atmosphere

37°C, 5% CO₂, 

Flask Coating


XXXXXX

Freezing Procedure

Shipping Conditions

Storage Conditions

—150 —196 °C



Sterility

PCR

