

4T1 | 300300

Tumorigenic BALB/c

Culture Medium RPMI 1640, w: 2.0 mM , w: 2.0 / NaHCO₃ (Cytion 820700a)

Supplements FBS 10%

Dissociation Reagent

Subculturing PBS T25

Freeze medium (FBS) + 10% DMSO

Thawing and Culturing Cells

1.
2. cryovial -150
3. 37°C
4. 70%
5. 15
6. 300 x g 3
7. 10
8.

Incubation Atmosphere 37°C, 5% CO₂,

Flask Coating

