

Product sheet

L Wnt-3A | 305184

Protein expression Wnt-3A

XXXXXXXXXX

Culture Medium DMEM, w: 4.5 g/l, w: 4 g/l, w: 3.7 g/l NaHCO₃, w: 1.0 g/l

Supplements FBS 10%, G-418 0.4 mg/ml

Dissociation Reagent XXXXXXXX

Subculturing PBS, T25

Fluid renewal 2-3

Freeze medium (FBS) + 10% DMSO

- Thawing and Culturing Cells**
1. XXXXXXXX
 2. XXXXXXXX cryovial -150 XXXXXXXX
 3. XXXXXXXX 37°C
 4. XXXXXXXX 70%
 5. XXXXXXXX 15
 6. XXXXXXXX 300 x g 3
 7. XXXXXXXX 10
 8. XXXXXXXX

Incubation Atmosphere 37°C, 5% CO₂, XXXXXXXX

