











Product sheet

Hep-CLS-E1  | 400196




Description	 C3H/HE  C3H/HE 
Organism	
Tissue	
Disease	
Synonyms	HEP-CLS-E1, E1




Breed/Subspecies	C3H/HE
Age	
Gender	
Morphology	
Growth properties	



















Citation	Hep-CLS-E1 ( Cytion 400196)
Biosafety level	1
NCBI_TaxID	10090
CellosaurusAccession	CVCL_5777



Protein expression	
---------------------------	---

Product sheet

Hep-CLS-E1  | 400196

Tumorigenic	  C3H/HE
Mutational profile	 53 
	
Culture Medium	Ham's F12, w: 1.0 mM  , w: 1.0 mM  , w: 1.1 g/L NaHCO ₃ ( Cytion 820600a)
Supplements	 FBS 10%
Dissociation Reagent	
Subculturing	 PBS   T25 
Seeding density	1×10^4  /  ^
Fluid renewal	 3  5 
Post-Thaw Recovery	 4  24 
Freeze medium	 ( FBS) + 10% DMSO

Hep-CLS-E1 | 400196

Thawing and Culturing Cells

- 1.
2. cryovial -150
3. 37°C
4. 70%
5. 15
6. 300 x g 3
7. 10
- 8.

Incubation Atmosphere

37°C, 5% CO₂,

Flask Coating

Freezing Procedure

Shipping Conditions

Storage Conditions

-150 -196 °C

/ / HLA

Sterility

PCR