

Product sheet

HB-CLS-2 [REDACTED] | 300191

[REDACTED]

<b>Description</b>	[REDACTED] ( [REDACTED] III) [REDACTED]
<b>Organism</b>	[REDACTED]
<b>Tissue</b>	[REDACTED]
<b>Disease</b>	[REDACTED]
<b>Synonyms</b>	HBCLS2

[REDACTED]

<b>Age</b>	50 [REDACTED]
<b>Gender</b>	[REDACTED]
<b>Ethnicity</b>	[REDACTED]
<b>Morphology</b>	[REDACTED]
<b>Growth properties</b>	[REDACTED]

[REDACTED]

<b>Citation</b>	HB-CLS-2 ([REDACTED] Cytion 300191)
<b>Biosafety level</b>	1
<b>NCBI_TaxID</b>	9606
<b>CellosaurusAccession</b>	CVCL_6214

[REDACTED]

<b>Tumorigenic</b>	[REDACTED]
--------------------	------------

[REDACTED]

**Product sheet**

**HB-CLS-2** | 300191

**Culture Medium** DMEM:Ham's F12 (1:1), w: 3.1 g/L, w: 2.5 g/L L-Asparagine, w: 15 g/L HEPES, w: 0.5 g/L

**Supplements** FBS 10%

**Dissociation Reagent**

**Subculturing** PBS, T25

**Seeding density**  $1 \times 10^4$  cells/cm<sup>2</sup>

**Fluid renewal** 2-3 times

**Post-Thaw Recovery** 24 hours

**Freeze medium** (FBS) + 10% DMSO

- Thawing and Culturing Cells**
1. Thaw cryovial in 37°C water bath
  2. Transfer cells to cryovial
  3. Centrifuge at 300 x g for 3 min
  4. Wash cells with PBS
  5. Resuspend cells in 15 µl PBS
  6. Centrifuge at 300 x g for 3 min
  7. Wash cells with PBS
  8. Resuspend cells in 10 µl PBS

**Incubation Atmosphere** 37°C, 5% CO<sub>2</sub>

