



CC531 | 500387

**Supplements** FBS 10% HEPES 20 mM

**Dissociation Reagent**

**Subculturing** PBS T25

**Seeding density** 1 2 x 10<sup>4</sup> 3 4

**Fluid renewal** 2 3

**Post-Thaw Recovery** 48

**Freeze medium** (FBS) + 10% DMSO

**Thawing and Culturing Cells**

- 1.
2. cryovial -150
3. 37°C
4. 70%
5. 15
6. 300 x g 3
7. 10
- 8.

**Incubation Atmosphere** 37°C, 5% CO<sub>2</sub>

**Flask Coating**

