

769-P | 300106

Thawing and Culturing Cells

- 1.
2. cryovial -150
3. 37°C
4. 70%
5. 15
6. 300 x g 3
7. 10
- 8.

Incubation Atmosphere 37°C, 5%_{CO2}

Flask Coating

Freezing Procedure

Shipping Conditions

Storage Conditions -150 -196 °C

/ / HLA

Sterility PCR

769-P ~~XXXXXXXX~~ | 300106

~~XXXXXXXXXXXX~~ HLA **A***: '03:01:01, '24:02:01
 B*: '07:02:01
 C*: '07:02:01
 DRB1*: 15:01:01G
 DQA1*: '01:02:01
 DQB1*: '06:02:01
 DPB1*: '04:01:01
 E: '01:03:02