

769-P [REDACTED] | 300106

Thawing and Culturing Cells

1. [REDACTED]
2. [REDACTED] cryovial [REDACTED]-150 [REDACTED]
3. [REDACTED] 37°C [REDACTED]
4. [REDACTED] 70%
5. [REDACTED] 15 [REDACTED]
6. [REDACTED] 300 x g [REDACTED] 3 [REDACTED]
7. [REDACTED] 10 [REDACTED]
8. [REDACTED]

Incubation Atmosphere 37°C, 5% CO₂, [REDACTED]

Flask Coating [REDACTED]

Freezing Procedure [REDACTED]


Shipping Conditions [REDACTED]

Storage Conditions [REDACTED]-150 [REDACTED]-196 °C [REDACTED]

[REDACTED] / [REDACTED] / HLA

Sterility [REDACTED] PCR [REDACTED]
[REDACTED]

769-P  | 300106

 HLA **A***: '03:01:01, '24:02:01
 B*: '07:02:01
 C*: '07:02:01
 DRB1*: 15:01:01G
 DQA1*: '01:02:01
 DQB1*: '06:02:01
 DPB1*: '04:01:01
 E: '01:03:02