

XXXXXXXXXXXXXXXXXXXX

XXXXXXXXXX

XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX













[illegible]

XXXXXX HS-729 | 300443

Isoenzymes	6 6, 6
Culture Medium	DMEM, w: 4.5 / 4.5, w: 4 4-4, w: 3.7 / 4.5 NaHCO3, w: 1.0 4
Supplements	4 FBS 10%
Dissociation Reagent	4
Subculturing	4 PBS 4 4
Seeding density	1×10^4 2/4
Fluid renewal	2 3 4
Post-Thaw Recovery	4 4 24 4
Freeze medium	4 (4 FBS) + 1

XXXXXX HS-729 | 300443

Thawing and Culturing Cells

1.  
2.   cryovial  -150  
3.   37°C 
4.   70%
5.   15 
6.  300 x g  3   
7.   10  
8. 

Incubation Atmosphere

37°C, 5% CO₂, 

Flask Coating

XXXXXX

Freezing Procedure

Shipping Conditions

Storage Conditions

-150 -196 °C

Sterility

PCR