





















MV4-11  | 300295

### Thawing and Culturing Cells

1. 
2.  cryovial  -150 
3.  37°C 
4.  70%
5.  15 
6.  300 x g  3 
7.  10 
8. 

### Incubation Atmosphere

37°C, 5% CO<sub>2</sub>, 

### Flask Coating



### Freezing Procedure



### Shipping Conditions




### Storage Conditions

 -150  -196 °C 

 /  / HLA

### Sterility

 PCR 

MV4-11 ~~XXXXXXXX~~ | 300295

---

~~XXXXXXXXXXXX~~ HLA    **A\***: '03:01:01, '68:01:02  
                                  **B\***: 14:02:01, 18:01:01  
                                  **C\***: '08:02:01, '15:02:01  
                                  **DRB1\***: '01:01:01, '13:02:01  
                                  **DQA1\***: '01:01:01, '01:02:01  
                                  **DQB1\***: '05:01:01, '06:09:01  
                                  **DPB1\***: '02:01:02, '04:01:01  
                                  **E**: '01:01, '01:03