



MV4-11  | 300295

Thawing and Culturing Cells

1. 
2.  cryovial  -150 
3.  37°C 
4.  70%
5.  15 
6.  300 x g  3 
7.  10 
8. 

Incubation Atmosphere

37°C, 5% CO₂, 

Flask Coating




Freezing Procedure



Shipping Conditions



Storage Conditions

 -150  -196 °C 

 /  / HLA

Sterility

 PCR 



MV4-11 ~~XXXXXXXX~~ | 300295

~~XXXXXXXXXXXX~~ HLA **A***: '03:01:01, '68:01:02
 B*: 14:02:01, 18:01:01
 C*: '08:02:01, '15:02:01
 DRB1*: '01:01:01, '13:02:01
 DQA1*: '01:01:01, '01:02:01
 DQB1*: '05:01:01, '06:09:01
 DPB1*: '02:01:02, '04:01:01
 E: '01:01, '01:03