







**XXXXXXXX Capan-1 | 300143**

---

**XXXXXXXXXXXX HLA**    **A\***: '01:01:01, '30:01:01  
                                  **B\***: 13:02:01, 57:01:01  
                                  **C\***: '06:02:01  
                                  **DRB1\***: '07:01:01, '13:05:01  
                                  **DQA1\***: '02:01:01, '05:05:01  
                                  **DQB1\***: '02:02:01, '03:01:01  
                                  **DPB1\***: '03:01:01G, '04:01:01G  
                                  **E**: '01:01:01