





**CaSki | 300145**

**Thawing and Culturing Cells**

1. [REDACTED]
2. [REDACTED] cryovial [REDACTED] -150 [REDACTED]
3. [REDACTED] 37°C [REDACTED]
4. [REDACTED] 70% [REDACTED]
5. [REDACTED] 15 [REDACTED]
6. [REDACTED] 300 x g [REDACTED] 3 [REDACTED]
7. [REDACTED] 10 [REDACTED]
8. [REDACTED]

**Incubation Atmosphere** 37°C, 5% CO<sub>2</sub>, [REDACTED]

**Flask Coating** [REDACTED]

**Freezing Procedure** [REDACTED]

**Shipping Conditions** [REDACTED]

**Storage Conditions** [REDACTED] -150 [REDACTED] -196 °C [REDACTED]

[REDACTED] / [REDACTED] / HLA

**Sterility** [REDACTED] PCR [REDACTED]  
[REDACTED]

XXXXXXXXXX CaSki | 300145

---

XXXXXXXXXXXX HLA    **A\***: '02:01:01, '03:01:01  
                          **B\***: '07:02:01, '37:01:01  
                          **C\***: '07:02:01  
                          **DRB1\***: '08:01:01G, '15:01:01G  
                          **DQA1\***: '01:02:01, '04:02  
                          **DQB1\***: '04:02:01, '06:02:01  
                          **DPB1\***: '04:01:01  
                          **E**: '01:03:02