

RAJI | 300359

Description	RAJI R.J.V. Pulvertaft 1963 23.2 CD19
-------------	---

Organism	
----------	--

Tissue	
--------	--

Disease	
---------	--

Synonyms	, P1-, GM04671
----------	----------------

Age	11
-----	----

Gender	
--------	--

Ethnicity	,
-----------	---

Cell type	
-----------	--

Growth properties	
-------------------	--

Citation	RAJI (Cytion 300359)
----------	-----------------------

Biosafety level	1
-----------------	---

NCBI_TaxID	9606
------------	------

CellosaurusAccession	CVCL_0511
----------------------	-----------

Products	DMEM (Gibco), FBS (Gibco), Penicillin/Streptomycin (Gibco), Fetal Bovine Serum (Gibco), RPMI 1640 (Gibco), Sodium Bicarbonate (Gibco), Cytion 820700a (Gibco)
Culture Medium	RPMI 1640, w: 2.0 mM NaHCO_3 (Gibco), w: 2.0 mM NaHCO_3 (Gibco) Cytion 820700a (Gibco)
Supplements	FBS 10%
Subculturing	5
Freeze medium	DMEM (Gibco), FBS (Gibco), Penicillin/Streptomycin (Gibco) (FBS) + 10% DMSO
Thawing and Culturing Cells	1. Thaw cells in a 37°C water bath. 2. Add cells to a cryovial and store at -150 °C. 3. Add cells to a 37°C water bath. 4. Add cells to a 70% FBS. 5. Add cells to a 15 mL. 6. Add cells to a 300 x g 13 mL. 7. Add cells to a 10 mL. 8. Add cells to a 10 mL.
Incubation Atmosphere	37°C, 5% CO_2 , 5% CO_2
Flask Coating	Flask Coating
Freezing Procedure	Freezing Procedure

RAJI | 300359

Shipping
Conditions

Storage
Conditions

Sterility

STR

CSF1PO: 10,12
D13S317: 13
D16S539: 8,11
D5S818: 10,13
D7S820: 10
TH01: 6,7
TPOX: 8,13
vWA: 16,19
D3S1358: 15,16
D21S11: 28,31
D18S51: 17
Penta E: 5,13
Penta D: 3,2,9
D8S1179: 14,15
FGA: 19,27

HLA

A*: '03:01:01
B*: 15:10:01
C*: '03:04:02, '04:01:01
DRB1*: '03:01:01, '10:01:01
DQA1*: '01:05:01, '05:01:01
DQB1*: '02:01:01, '05:01:01
DPB1*: '01:01:01
E: '01:01:01