

L6565 | 305189

Description	L6565 L6565
-------------	----------------

Organism	
----------	--

Tissue	
--------	--

Morphology	
------------	--

Growth properties	
-------------------	--

Citation	L6565 (Cytion 305189)
----------	------------------------

Biosafety level	1
-----------------	---

NCBI_TaxID	10090
------------	-------

CellosaurusAccession	CVCL_A9NB
----------------------	-----------

Culture Medium	DMEM:Ham's F12 (1:1), w: 3.1 / , w: 2.5 L-, w: 15 HEPES, w: 0.5
----------------	---

Supplements	FBS 10%, 0.005 / , 0.01 / ,
-------------	-----------------------------

Subculturing	5
--------------	---

Fluid renewal	2 3
---------------	-----

L6565 | 305189

Freeze
medium

DMSC medium (100% FBS) + 10% DMSO

Thawing and
Culturing Cells

1. Thaw the cells in a water bath at 37°C.
2. Add the cells to a cryovial and store at -150 °C.
3. Thaw the cells in a water bath at 37°C.
4. Add the cells to a flask and incubate at 37°C.
5. Add the cells to a flask and incubate at 37°C.
6. Add the cells to a flask and incubate at 37°C.
7. Add the cells to a flask and incubate at 37°C.
8. Add the cells to a flask and incubate at 37°C.

Incubation
Atmosphere

37°C, 5% CO₂, humidified

Flask Coating

None

Freezing
Procedure

None

Shipping
Conditions

None

Storage
Conditions

None

L6565 | 305189

Genomic map of the 100 kb region around the PCRC gene. The top track shows the PCRC gene structure with exons as boxes and introns as lines. The bottom track shows the genomic context with various features and annotations.