

XXXXXXXXXXXXXXXXXXXXXXXXXXXX

XXXXXXXXXX

(X)(X)(X)(X)(X)(X)(X)(X)(X)(X)(X)(X)(X)

XXXXXXXXXXXXXXXXXXXX

Product sheet

XXXXXTERA-C3H | 400360

Culture Medium	DMEM, w: 4.5 g/l / 100 ml, w: 4 g/l - 100 ml, w: 3.7 g/l / 100 ml NaHCO3, w: 1.0 g/l - 100 ml
Supplements	XXXXX FBS 10%
Dissociation Reagent	XXXXX
Subculturing	XXXXX PBS XXXXX T25 XXXX
Fluid renewal	2 XXXX
Freeze medium	XXXXX (XXXXX FBS) + 10% DMSO
Thawing and Culturing Cells	1. XXXX 2. XXXX cryovial XXXX -150 XXXX 3. XXXX 37°C XXXX 4. XXXX 70% 5. XXXX 15 XXXX 6. XXXX 300 x g XXXX 3 XXXX 7. XXXX 10 XXXX 8. XXXX
Incubation Atmosphere	37°C, 5% CO2 XXXX
Flask Coating	XXXXX
Freezing Procedure	XXXXX

XXXXXX TERA-C3H | 400360

Shipping Conditions

Storage Conditions

Sterility