

TM3 | 305167

**Description** TM3 Leydig 11 13

**Organism**

**Tissue**

**Synonyms** TM-3

**Breed/Subspecies** BALB/c

**Age** 11 13

**Gender**

**Morphology**

**Growth properties**

**Citation** TM3 (Cytion 305167)

**Biosafety level** 1


**NCBI\_TaxID** 10090

**CellosaurusAccession** CVCL\_4326



**Culture Medium** DMEM:Ham's F12 (1:1), w: 3.1 / w: 2.5 L- w: 15 HEPES, w: 0.5




TM3 | 305167

**Supplements**  FBS 2.5%  5%





















**Dissociation Reagent** 


**Subculturing**  PBS   T25 

**Fluid renewal** 2  3 

**Freeze medium**   (  FBS ) + 10% DMSO

**Thawing and Culturing Cells**

1.  
2.   cryovial  -150 
3.   37°C 
4.   70%
5.   15 
6.  300 x g  3 
7.  10 
8. 

**Incubation Atmosphere** 37°C, 5% CO<sub>2</sub>, 

**Flask Coating** 

**Freezing Procedure** 

TM3 | 305167

**Shipping Conditions** [REDACTED]

**Storage Conditions**

[REDACTED] -150 [REDACTED] -196 °C [REDACTED]

/ / HLA

**Sterility** [REDACTED] PCR [REDACTED]  
[REDACTED] [REDACTED] [REDACTED]