

WPMY-1 | 305083

Thawing and Culturing Cells

1.

...
2.

... cryovial ... -150 ...
3.

... 37°C ...
4.

... 70% ...
5.

... 15 ...
6.

... 300 x g ... 3 ...
7.

... 10 ...
8.

...

Incubation Atmosphere

37°C, 5% CO₂, ...

Flask Coating

...

Freezing Procedure

...

Shipping Conditions

...

Storage Conditions

... -150 ... -196 °C ...

... / ... / HLA

Sterility

... PCR ...

...