

BJAB | 302006

Thawing and Culturing Cells

1. [REDACTED]
2. [REDACTED] cryovial [REDACTED]-150 [REDACTED]
3. [REDACTED] 37°C [REDACTED]
4. [REDACTED] 70%
5. [REDACTED] 15 [REDACTED]
6. [REDACTED] 300 x g [REDACTED] 3 [REDACTED]
7. [REDACTED] 10 [REDACTED]
8. [REDACTED]

Incubation Atmosphere 37°C, 5% CO₂, [REDACTED]

Flask Coating [REDACTED]

Freezing Procedure [REDACTED]

Shipping Conditions [REDACTED]

Storage Conditions [REDACTED]-150 [REDACTED]-196 °C [REDACTED]

[REDACTED] / [REDACTED] / HLA

Sterility [REDACTED] PCR [REDACTED]
[REDACTED] [REDACTED] [REDACTED]

XXXXXXXXBJAB | 302006

XXXXXXXXXXXX HLA **A***: 01:01:83, '02:01:01
 B*: 13:02:01, 35:01:01
 C*: '04:01:01, '06:02:01
 DRB1*: '12:01:01, '13:02:01
 DQA1*: '01:02:01, '05:05:01
 DQB1*: '03:01, '06:04:01
 DPB1*: '04:02:01G
 E: '01:01, '01:03