

LM/TK(LMTK-) | 305176

XXXXXXXXXX

Culture Medium DMEM, w: 4.5 g/l / 45 ml, w: 4 g/l - , w: 3.7 g/l / 37 ml NaHCO₃, w: 1.0 g/l

Supplements FBS 10%

Dissociation Reagent XXXXXXXX

Subculturing PBS XXXXXXXX T25

Fluid renewal 2 XXXXXXXX

Freeze medium XXXXXXXX (FBS) + 10% DMSO

- Thawing and Culturing Cells**
1. XXXXXXXX
 2. XXXXXXXX cryovial -150 XXXXXXXX
 3. XXXXXXXX 37°C
 4. XXXXXXXX 70%
 5. XXXXXXXX 15
 6. XXXXXXXX 300 x g 3 XXXXXXXX
 7. XXXXXXXX 10 XXXXXXXX
 8. XXXXXXXX

Incubation Atmosphere 37°C, 5% CO₂, XXXXXXXX

Flask Coating XXXXXXXX

