

HL-60 | 300209

Thawing and Culturing Cells

1. [REDACTED]
2. [REDACTED] cryovial [REDACTED] -150 [REDACTED]
3. [REDACTED] 37°C [REDACTED]
4. [REDACTED] 70%
5. [REDACTED] 15 [REDACTED]
6. [REDACTED] 300 x g [REDACTED] 3 [REDACTED]
7. [REDACTED] 10 [REDACTED]
8. [REDACTED]

Incubation Atmosphere

37°C, 5% CO₂, [REDACTED]

Flask Coating

[REDACTED]

Freezing Procedure

[REDACTED]

Shipping Conditions

[REDACTED]

Storage Conditions

[REDACTED] -150 [REDACTED] -196 °C [REDACTED]

[REDACTED] / [REDACTED] / HLA

Sterility

[REDACTED] PCR [REDACTED]

[REDACTED] [REDACTED] [REDACTED]

██████████ HL-60 | 300209

██████████████████ HLA **A***: '01:01:01
 B*: '57:01:01
 C*: '06:02:01
 DRB1*: '07:01:01
 DQA1*: '02:01:01
 DQB1*: '03:03:02
 DPB1*: '04:01:01, '13:01:01
 E: '01:01:01, '01:09