

BRL | 305193

Dissociation Reagent [REDACTED]

Subculturing [REDACTED] PBS [REDACTED] [REDACTED] T25 [REDACTED]

Fluid renewal 2 [REDACTED] 3 [REDACTED]

Freeze medium [REDACTED] ([REDACTED] FBS) + 10% DMSO

- Thawing and Culturing Cells**
1. [REDACTED]
 2. [REDACTED] cryovial [REDACTED] -150 [REDACTED]
 3. [REDACTED] 37°C [REDACTED]
 4. [REDACTED] 70%
 5. [REDACTED] 15 [REDACTED]
 6. [REDACTED] 300 x g [REDACTED] 3 [REDACTED]
 7. [REDACTED] 10 [REDACTED]
 8. [REDACTED]

Incubation Atmosphere 37°C, 5% CO₂, [REDACTED]

Flask Coating [REDACTED]

Freezing Procedure [REDACTED]

Shipping Conditions [REDACTED]

