





T84 | 300354

**Thawing and Culturing Cells**

1. [REDACTED]
2. [REDACTED] cryovial [REDACTED] -150 [REDACTED]
3. [REDACTED] 37°C [REDACTED]
4. [REDACTED] 70% [REDACTED]
5. [REDACTED] 15 [REDACTED]
6. [REDACTED] 300 x g [REDACTED] 3 [REDACTED]
7. [REDACTED] 10 [REDACTED]
8. [REDACTED]

**Incubation Atmosphere**

37°C, 5% CO<sub>2</sub>, [REDACTED]

**Flask Coating**

[REDACTED]

**Freezing Procedure**

[REDACTED]

**Shipping Conditions**

[REDACTED]

**Storage Conditions**

[REDACTED] -150 [REDACTED] -196 °C [REDACTED]

/ / HLA

**Sterility**

[REDACTED] PCR [REDACTED]

[REDACTED]

T84 | 300354

---

**HLA**    **A\***: '02:01:01, '24:02:01  
          **B\***: '18:01:01, '35:01:01  
          **C\***: '04:01:01, '07:01:01  
          **DRB1\***: '01:01:01, '09:01:02  
          **DQA1\***: '01:01:01, '03:02:01  
          **DQB1\***: '03:03:02, '05:01:01  
          **DPB1\***: '02:01:02, '04:01:01  
          **E**: '01:03:01, '01:03:02