

RPMI 1640, w: 2.0 mM , w: 2.0 / N
aHCO3 | 820700a

RPMI 1640 Medium RPMI medium Roswell Park Memorial Institute (RPMI)

RPMI 1640 Medium

RPMI 1640 Medium Serum FBS) 10%

RPMI 1640 Medium CO2 5-10% pH

-
- +2°C +8°C
- 4°C 6-8
-
- +2°C +8°C +37°C
- 37°C
-

	()
	10.00
	434.40
	200.00
	56.82
	20.00
	65.20
	20.00
	20.27

, w: 2.0 / N

	001- 00000000- 000- 000000	20.00
	001- 0000000000	50.00
	001- 000000	50.00
	001- 00000 000000000000000000	40.00
	001- 00000000	15.00
	001- 000000000000	15.00
	001- 000000	20.00
	001- 00000	30.00
	001- 0000000	20.00
	001- 00000000	5.00
	001- 00000000 2Na 2H2O	28.83
	001- 00000	20.00
000000	0000000- 000000000000	1.00
	00- 000000	0.20
	0000000000	3.00
	00- 000000000000000000	0.25
	0000000	1.00
	0000- 00000000	35.00
	0000000000	1.00
	00000000 0000000000000000	1.00
	000000000	0.20

Product sheet

RPMI 1640, w: 2.0 mM, w: 2.0 / N
aHCO3 | 820700a

	Glucose	1.00
	Insulin	0.005
Transferrin	Transferrin	100.00
	Sealant	400.00
	Penicillin	100.00
	NaCl	6000.00
	Streptomycin	2000.00
	Na2HPO4	800.00
Hydrocortisone	Hydrocortisone	2000.00
	Hydrocortisone	1.00
	Hydrocortisone	5.30