





**XXXXXXXXGCT | 300155**

**Flask Coating**

**Freezing Procedure**

**Shipping Conditions**

**Storage Conditions**

-150  -196 °C

**XXXXXXXXXXXXXXXXXXXX / XXXXXXXXXXXXXXXXXXXXXXX / HLA**

**Sterility**  PCR

- XXXXXXXXXXXX HLA**
- A\*:** '01:01:01, '23:01:01
  - B\*:** '08:01:01, '15:17:01
  - C\*:** '07:01:01, '07:01:02
  - DRB1\*:** '03:01:01, '04:04:01
  - DQA1\*:** '03:01:01, '05:01:01
  - DQB1\*:** '02:01:01, '03:02:01
  - DPB1\*:** '01:01:01, '02:01:02
  - E:** '01:01:01, '01:03:05