





3T3-L1 | 400107

Thawing and Culturing Cells

1. [REDACTED]
2. [REDACTED] cryovial [REDACTED]-150 [REDACTED]
3. [REDACTED] 37°C [REDACTED]
4. [REDACTED] 70%
5. [REDACTED] 15 [REDACTED]
6. [REDACTED] 300 x g [REDACTED] 3 [REDACTED]
7. [REDACTED] 10 [REDACTED]
8. [REDACTED]

**Incubation Atmosphere** 37°C, 5% CO<sub>2</sub>, [REDACTED]

**Flask Coating** [REDACTED]

**Freezing Procedure** [REDACTED]

**Shipping Conditions** [REDACTED]

**Storage Conditions** [REDACTED]-150 [REDACTED]-196 °C [REDACTED]

[REDACTED] / [REDACTED] / HLA

**Sterility** [REDACTED] PCR [REDACTED]  
 [REDACTED] [REDACTED] [REDACTED]