

A498 | 300113

Thawing and Culturing Cells

- 1.
2. cryovial -150
3. 37°C
4. 70%
5. 15
6. 300 x g 3
7. 10
- 8.

Incubation Atmosphere 37°C, 5%_{CO2},

Flask Coating

Freezing Procedure

Shipping Conditions

Storage Conditions -150 -196 °C

/ / HLA

Sterility PCR

A498  | 300113

 HLA **A***: '02:01:01
 B*: '08:01:01
 C*: '07:01:01
 DRB1*: '03:01:01
 DQA1*: '05:01:01
 DQB1*: '02:01:01
 DPB1*: '01:01:01
 E: '01:03:02