

Product sheet

D341Med | 305136

Protein expression [REDACTED], [REDACTED], [REDACTED], [REDACTED]

Tumorigenic [REDACTED]

Culture Medium EMEM (MEM Eagle), w: 2 mM L-Glutamine, w: 2.2 g/L NaHCO₃, w: EBSS ([REDACTED] Cytion 820100a)

Supplements [REDACTED] FBS 10% [REDACTED] NEAA 1%

Doubling time 37 [REDACTED]

Subculturing [REDACTED] 5 [REDACTED]

Fluid renewal 2 [REDACTED] 3 [REDACTED]

Freeze medium [REDACTED] [REDACTED] ([REDACTED] FBS) + 10% DMSO

- Thawing and Culturing Cells**
1. [REDACTED]
 2. [REDACTED] cryovial [REDACTED] -150 [REDACTED] [REDACTED]
 3. [REDACTED] [REDACTED] 37°C [REDACTED]
 4. [REDACTED] [REDACTED] 70%
 5. [REDACTED] [REDACTED] 15 [REDACTED]
 6. [REDACTED] 300 x g [REDACTED] 3 [REDACTED] [REDACTED]
 7. [REDACTED] [REDACTED] 10 [REDACTED] [REDACTED]
 8. [REDACTED]

Incubation Atmosphere 37°C, 5% CO₂, [REDACTED]

Product sheet

D341Med | 305136

Flask Coating [REDACTED]

Freezing Procedure

[REDACTED]

Shipping Conditions [REDACTED]

Storage Conditions

[REDACTED] -150 [REDACTED] -196 °C [REDACTED]

/ / HLA

Sterility [REDACTED] PCR [REDACTED]
[REDACTED] [REDACTED] [REDACTED]