



Product sheet

DSL-6A-C1  | 500166



**Tumorigenic**  Lewis 



**Culture Medium** Waymouth medium ()

**Supplements**  FBS 10%, L-glutamine 2.0 mM

**Dissociation Reagent** 

**Subculturing**  PBS  T25 

**Seeding density**  $1 \times 10^4$  

**Fluid renewal** 2 

**Post-Thaw Recovery**  24 

**Freeze medium**  (  FBS) + 10% DMSO

DSL-6A-C1 | 500166

### Thawing and Culturing Cells

- 1.
2. cryovial -150
3. 37°C
4. 70%
5. 15
6. 300 x g 3
7. 10
- 8.

### Incubation Atmosphere

37°C, 5% CO<sub>2</sub>,

### Flask Coating

### Freezing Procedure

### Shipping Conditions

### Storage Conditions

-150 -196 °C

/ / HLA

### Sterility

PCR