

K7M2 wt | 305188

Receptors expressed (C3), Fc, IgG, I(Fcgr1),

Tumorigenic

Culture Medium DMEM, w: 4.5 / w: 4 - w: 3.7 / NaHCO₃, w: 1.0

Supplements FBS 10%

Dissociation Reagent

Subculturing PBS T25

Fluid renewal 2 3

Freeze medium (FBS) + 10% DMSO

- Thawing and Culturing Cells**
1.
 2. cryovial -150
 3. 37°C
 4. 70%
 5. 15
 6. 300 x g 3
 7. 10
 8.

Product sheet

K7M2 wt | 305188

Incubation Atmosphere 37°C, 5% CO₂, [REDACTED]

Flask Coating [REDACTED]

Freezing Procedure [REDACTED]

Shipping Conditions [REDACTED]

Storage Conditions [REDACTED] -150 [REDACTED] -196 °C [REDACTED]

/ / HLA

Sterility [REDACTED] PCR [REDACTED]
[REDACTED] [REDACTED] [REDACTED]